

Product Datasheet

Human Xin Actin-binding Repeat-containing Protein 1 ELISA Kit (orb1209887)

Description	Human Xin Actin-binding Repeat-containing Protein 1, XIRP1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	XIRP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	198,561
Uniprot ID	Q702N8
Solubility (25°C)	4800ng/L
Sensitivity	9.78ng/L
Entrez	165904
Expiration Date	Please enquire.