

Product Datasheet

Human Fructose-bisphosphate Aldolase C ELISA Kit (orb1209900)

Description	Human Fructose-bisphosphate Aldolase C, ALDOC ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	ALDOC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,456
Uniprot ID	P09972
Solubility (25°C)	4800ng/L
Sensitivity	9.87ng/L
Entrez	230
Expiration Date	Please enquire.