

Product Datasheet

Human Complement C1Q-like Protein 3 ELISA Kit (orb1209906)

Description	Human Complement C1Q-like Protein 3, C1QL3 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	C1QL3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,719
Uniprot ID	Q5VWW1
Solubility (25°C)	24ng/ml
Sensitivity	0.022ng/ml
Entrez	389941
Expiration Date	Please enquire.