

Product Datasheet

Human Arf-Gap with Sh3 Domain ELISA Kit (orb1209928)

Description	Human Arf-Gap with Sh3 Domain, Ank Repeat and Ph Domain- Containing Protein 1, ASAP1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	ASAP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	125,498
Uniprot ID	Q9ULH1
Solubility (25°C)	4800ng/L
Sensitivity	9.91ng/L
Entrez	50807
Expiration Date	Please enquire.