

Product Datasheet

Human Phosphatidylinositol-glycan-specific Phospholipase D ELISA Kit (orb1209936)

Description	Human Phosphatidylinositol-glycan-specific Phospholipase D, GPLD1 ELISA Kit
Reactivity	Human
Range	1-400ng/mL
Tested Applications	ELISA
Target	GPLD1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	92,336
Uniprot ID	P80108
Solubility (25°C)	480ng/mL
Sensitivity	0.53ng/mL
Entrez	2822
Expiration Date	Please enquire.