

Product Datasheet

Human Uncharacterized Protein C9Orf72 ELISA Kit (orb1209948)

Description	Human Uncharacterized Protein C9Orf72, C9ORF72 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	C9orf72
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,481
Uniprot ID	A6NCS6
Solubility (25°C)	4800ng/L
Sensitivity	10.14ng/L
Entrez	257407
Expiration Date	Please enquire.