

Product Datasheet

Human At Rich Interactive Domain Containing Protein 1A ELISA Kit (orb1209952)

Description	Human At Rich Interactive Domain Containing Protein 1A, ARID1A ELISA Kit
Reactivity	Human
Range	12.5-800ng/L
Tested Applications	ELISA
Target	ARID1A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	242,045
Uniprot ID	O14497
Solubility (25°C)	1600ng/L
Sensitivity	6.61ng/L
Entrez	8289
Expiration Date	Please enquire.