

Product Datasheet

Human Ubiquitin Carboxyl Terminalhydrolase Bap1 ELISA Kit (orb1209954)

Description	Human Ubiquitin Carboxyl Terminalhydrolase Bap1, BAP1 ELISA Kit
Reactivity	Human
Range	20-7000ng/ml
Tested Applications	ELISA
Target	BAP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	80,362
Uniprot ID	Q92560
Solubility (25°C)	8000ng/ml
Sensitivity	9.91ng/ml
Entrez	8314
Expiration Date	Please enquire.