

Product Datasheet

Human Serine Threonine-protein Kinase B-raf ELISA Kit (orb1209958)

Description	Human Serine Threonine-protein Kinase B-raf, BRAF ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	BRAF
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	84,437
Uniprot ID	P15056
Solubility (25°C)	4800ng/L
Sensitivity	9.11ng/L
Entrez	673
Expiration Date	Please enquire.