

Product Datasheet

Human Low Affinity Immunoglobulin Gamma Fc Region Receptor 2B ELISA Kit (orb1209959)

Description	Human Low Affinity Immunoglobulin Gamma Fc Region Receptor 2B, FCGR2B ELISA Kit
Reactivity	Human
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	FCGR2B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	34,044
Uniprot ID	P31994
Solubility (25°C)	96ng/ml
Sensitivity	0.13ng/ml
Entrez	2213
Expiration Date	Please enquire.