

Product Datasheet

Human Ferrochelatase ELISA Kit (orb1209964)

Description	Human Ferrochelatase, Mitochondrial, FECH ELISA Kit
Reactivity	Human
Range	0.05-35ng/ml
Tested Applications	ELISA
Target	FECH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	47,862
Uniprot ID	P22830
Solubility (25°C)	40ng/ml
Sensitivity	0.028ng/ml
Entrez	2235
Expiration Date	Please enquire.