

Product Datasheet

Human Calcyphosine ELISA Kit (orb1209967)

Description	Human Calcyphosine, CAPS ELISA Kit
Reactivity	Human
Range	0.05-30mmol/L
Tested Applications	ELISA
Target	CAPS
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,240
Uniprot ID	Q13938
Solubility (25°C)	40mmol/L
Sensitivity	0.026mmol/L
Entrez	828
Expiration Date	Please enquire.