

Product Datasheet

Human Succinate Dehydrogenase Iron-sulfur Subunit ELISA Kit (orb1209972)

Description	Human Succinate Dehydrogenase Iron-sulfur Subunit, Mitochondrial, SDHB ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	SDHB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	31,630
Uniprot ID	P21912
Solubility (25°C)	40ng/ml
Sensitivity	0.027ng/ml
Entrez	6390
Expiration Date	Please enquire.