

Product Datasheet

Human Ectonucleotide Pyrophosphatase ELISA Kit (orb1209973)

Description	Human Ectonucleotide Pyrophosphatase, Phosphodiesterase Family Member 1, ENPP1 ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	ENPP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	104,924
Uniprot ID	P22413
Solubility (25°C)	12.8ng/ml
Sensitivity	0.023ng/ml
Entrez	5167
Expiration Date	Please enquire.