

Product Datasheet

Human Basic Salivary Proline Rich Protein 1 ELISA Kit (orb1209979)

Description	Human Basic Salivary Proline Rich Protein 1, PRB1 ELISA Kit
Reactivity	Human
Range	0.2-70ng/ml
Tested Applications	ELISA
Target	PRB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	38,562
Uniprot ID	P04280
Solubility (25°C)	80ng/ml
Sensitivity	0.12ng/ml
Entrez	5542
Expiration Date	Please enquire.