

Product Datasheet

Human Atp Synthase Subunit Beta ELISA Kit (orb1209981)

Description	Human Atp Synthase Subunit Beta, Mitochondrial, ATP5B ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	ATP5B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	56,560
Uniprot ID	P06576
Solubility (25°C)	3200ng/L
Sensitivity	7.32ng/L
Entrez	506
Expiration Date	Please enquire.