

Product Datasheet

Human Acylamino-acid-releasing Enzyme ELISA Kit (orb1209982)

Description	Human Acylamino-acid-releasing Enzyme, APEH ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	APEH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	81,225
Uniprot ID	P13798
Solubility (25°C)	4800ng/L
Sensitivity	9.88ng/L
Entrez	327
Expiration Date	Please enquire.