

Product Datasheet

Human Fibroblast Growth Factor-20 ELISA Kit (orb1209989)

Description	Human Fibroblast Growth Factor-20, FGF20 ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	FGF-20
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,499
Uniprot ID	Q9NP95
Solubility (25°C)	1600ng/L
Sensitivity	3.98ng/L
Entrez	26281
Expiration Date	Please enquire.