

Product Datasheet

Human A Disintegrin and Metalloproteinase with Thrombospondin Motifs 8 ELISA Kit (orb1209992)

Description	Human A Disintegrin and Metalloproteinase with Thrombospondin Motifs 8, ADAMTS8 ELISA Kit
Reactivity	Human
Range	0.5-200ng/mL
Tested Applications	ELISA
Target	ADAMTS8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	96,460
Uniprot ID	Q9UP79
Solubility (25°C)	240ng/mL
Sensitivity	0.27ng/mL
Entrez	11095
Expiration Date	Please enquire.