

Product Datasheet

Human Spondin-2 ELISA Kit (orb1210003)

Description	Human Spondin-2, SPON2 ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	Spon2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,846
Uniprot ID	Q9BUD6
Solubility (25°C)	40ng/ml
Sensitivity	0.031ng/ml
Entrez	10417
Expiration Date	Please enquire.