

Product Datasheet

Human Interferon Regulatory Factor 7 ELISA Kit (orb1210014)

Description	Human Interferon Regulatory Factor 7, IRF7 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	IRF7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	54,278
Uniprot ID	Q92985
Solubility (25°C)	3200ng/L
Sensitivity	7.76ng/L
Entrez	3665
Expiration Date	Please enquire.