

Product Datasheet

Human Interferon Regulatory Factor 6 ELISA Kit (orb1210015)

Description	Human Interferon Regulatory Factor 6, IRF6 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	IRF6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	53,130
Uniprot ID	014896
Solubility (25°C)	2400ng/L
Sensitivity	5.78ng/L
Entrez	3664
Expiration Date	Please enquire.