

Product Datasheet

Human T-cell Surface Glycoprotein Cd3 Zeta Chain ELISA Kit (orb1210018)

Description	Human T-cell Surface Glycoprotein Cd3 Zeta Chain, CD247 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	CD247
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,696
Uniprot ID	P20963
Solubility (25°C)	24ng/ml
Sensitivity	0.026ng/ml
Entrez	919
Expiration Date	Please enquire.