

Product Datasheet

Human Aldo-Keto Reductase Family 1 Member C3 ELISA Kit (orb1210026)

Description	Human Aldo-Keto Reductase Family 1 Member C3, AKR1C3 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	AKR1C3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,853
Uniprot ID	P42330
Solubility (25°C)	3200ng/L
Sensitivity	8.32ng/L
Entrez	8644
Expiration Date	Please enquire.