

## Product Datasheet

# Human Neurogenic Locus Notch Homolog Protein 3 ELISA Kit (orb1210029)

|                            |                                                                  |
|----------------------------|------------------------------------------------------------------|
| <b>Description</b>         | Human Neurogenic Locus Notch Homolog Protein 3, NOTCH3 ELISA Kit |
| <b>Reactivity</b>          | Human                                                            |
| <b>Range</b>               | 70-15000ng/L                                                     |
| <b>Tested Applications</b> | ELISA                                                            |
| <b>Target</b>              | NOTCH3                                                           |
| <b>Storage</b>             | 2-8°C                                                            |
| <b>Note</b>                | For research use only                                            |
| <b>Application notes</b>   | Sandwich                                                         |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates                           |
| <b>MW</b>                  | 243,631                                                          |
| <b>Uniprot ID</b>          | <b>Q9UM47</b>                                                    |
| <b>Solubility (25°C)</b>   | 16000ng/L                                                        |
| <b>Sensitivity</b>         | 36.17ng/L                                                        |
| <b>Entrez</b>              | <b>4854</b>                                                      |
| <b>Expiration Date</b>     | Please enquire.                                                  |