

Product Datasheet

Human Deoxyribonuclease-1-like 1 ELISA Kit (orb1210038)

Description	Human Deoxyribonuclease-1-like 1, DNASE1L1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	DNASE1L1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	33,893
Uniprot ID	P49184
Solubility (25°C)	3200ng/L
Sensitivity	7.77ng/L
Entrez	1774
Expiration Date	Please enquire.