

Product Datasheet

Human Endothelin-converting Enzyme 1 ELISA Kit (orb1210039)

Description	Human Endothelin-converting Enzyme 1, ECE1 ELISA Kit
Reactivity	Human
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	ECE1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	87,164
Uniprot ID	P42892
Solubility (25°C)	96ng/ml
Sensitivity	0.14ng/ml
Entrez	1889
Expiration Date	Please enquire.