

Product Datasheet

Human Tubulin Polymerization-promoting Protein Family Member 3 ELISA Kit (orb1210042)

Description	Human Tubulin Polymerization-promoting Protein Family Member 3, TPPP3 ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	TPPP3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,985
Uniprot ID	Q9BW30
Solubility (25°C)	12.8ng/ml
Sensitivity	0.017ng/ml
Entrez	51673
Expiration Date	Please enquire.