

Product Datasheet

Human Spry Domain-containing Socs Box Protein 2 ELISA Kit (orb1210048)

Description	Human Spry Domain-containing Socs Box Protein 2, SPSB2 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	SPSB2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,630
Uniprot ID	Q99619
Solubility (25°C)	2400ng/L
Sensitivity	5.12ng/L
Entrez	84727
Expiration Date	Please enquire.