

Product Datasheet

Human Isocitrate Dehydrogenase ELISA Kit (orb1210054)

Description	Human Isocitrate Dehydrogenase, Nadp Cytoplasmic, IDH1 ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	IDH1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,659
Uniprot ID	075874
Solubility (25°C)	80ng/ml
Sensitivity	0.11ng/ml
Entrez	3417
Expiration Date	Please enquire.