

Product Datasheet

Human Follistatin-related Protein 5 ELISA Kit (orb1210060)

Description	Human Follistatin-related Protein 5, FSTL5 ELISA Kit
Reactivity	Human
Range	0.5-150ng/ml
Tested Applications	ELISA
Target	FSTL5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	95,751
Uniprot ID	Q8N475
Solubility (25°C)	160ng/ml
Sensitivity	0.27ng/ml
Entrez	56884
Expiration Date	Please enquire.