

Product Datasheet

Human Centromere Protein H ELISA Kit (orb1210067)

Description	Human Centromere Protein H, CENPH ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	CENPH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,481
Uniprot ID	Q9H3R5
Solubility (25°C)	48ng/ml
Sensitivity	0.054ng/ml
Entrez	64946
Expiration Date	Please enquire.