

Product Datasheet

Human Surfeit Locus Protein 2 ELISA Kit (orb1210070)

Description	Human Surfeit Locus Protein 2, SURF2 ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	SURF2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,648
Uniprot ID	Q15527
Solubility (25°C)	4000ng/L
Sensitivity	9.44ng/L
Entrez	6835
Expiration Date	Please enquire.