

Product Datasheet

Human D-amino Acidoxidase Activator ELISA Kit (orb1210072)

Description	Human D-amino Acidoxidase Activator, DAOA ELISA Kit
Reactivity	Human
Range	0.1-30ng/ml
Tested Applications	ELISA
Target	DAOA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,108
Uniprot ID	P59103
Solubility (25°C)	32ng/ml
Sensitivity	0.051ng/ml
Entrez	267012
Expiration Date	Please enquire.