

Product Datasheet

Human Microtubule Associated Protein RP/EB Family Member 1 ELISA Kit (orb1210076)

Description	Human Microtubule Associated Protein RP/EB Family Member 1, MAPRE1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	MAPRE1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,999
Uniprot ID	Q15691
Solubility (25°C)	3200ng/L
Sensitivity	7.71ng/L
Entrez	22919
Expiration Date	Please enquire.