

Product Datasheet

Human Cochlin ELISA Kit (orb1210088)

Description	Human Cochlin, COCH ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	COCH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	59,483
Uniprot ID	O43405
Solubility (25°C)	24ng/ml
Sensitivity	0.03ng/ml
Entrez	1690
Expiration Date	Please enquire.