

Product Datasheet

Human Chromobox Protein Homolog 3 ELISA Kit (orb1210089)

Description	Human Chromobox Protein Homolog 3, CBX3 ELISA Kit
Reactivity	Human
Range	5-1000ng/L
Tested Applications	ELISA
Target	CBX3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,811
Uniprot ID	Q13185
Solubility (25°C)	1280ng/L
Sensitivity	2.22ng/L
Entrez	11335
Expiration Date	Please enquire.