

Product Datasheet

Human Carbohydrate Sulfotransferase 7 ELISA Kit (orb1210097)

Description	Human Carbohydrate Sulfotransferase 7, CHST7 ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	CHST7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	54,266
Uniprot ID	Q9NS84
Solubility (25°C)	16ng/ml
Sensitivity	0.021ng/ml
Entrez	56548
Expiration Date	Please enquire.