

Product Datasheet

Human Pituitary Adenylate Cyclase-activating Polypeptide ELISA Kit (orb1210120)

Description	Human Pituitary Adenylate Cyclase-activating Polypeptide, ADCYAP1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	ADCYAP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,835
Uniprot ID	P18509
Solubility (25°C)	3200ng/L
Sensitivity	6.87ng/L
Entrez	116
Expiration Date	Please enquire.