

Product Datasheet

Human Cystatin-F ELISA Kit (orb1210121)

Description	Human Cystatin-F, CST7 ELISA Kit
Reactivity	Human
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	CST7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,454
Uniprot ID	076096
Solubility (25°C)	96ng/ml
Sensitivity	0.12ng/ml
Entrez	8530
Expiration Date	Please enquire.