

Product Datasheet

Human Cystatin-S ELISA Kit (orb1210123)

Description	Human Cystatin-S, CST4 ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	CST4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,214
Uniprot ID	P01036
Solubility (25°C)	16ng/ml
Sensitivity	0.015ng/ml
Entrez	1472
Expiration Date	Please enquire.