

Product Datasheet

Human Lipocalin 15 ELISA Kit (orb1210131)

Description	Human Lipocalin 15, LCN15 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	LCN15
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,454
Uniprot ID	Q6UWW0
Solubility (25°C)	2400ng/L
Sensitivity	4.66ng/L
Entrez	389812
Expiration Date	Please enquire.