

Product Datasheet

Human Intellectin 2 ELISA Kit (orb1210145)

Description	Human Intellectin 2, ITLN2 ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	ITLN2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,212
Uniprot ID	Q8WWU7
Solubility (25°C)	6400ng/L
Sensitivity	10.13ng/L
Entrez	142683
Expiration Date	Please enquire.