

Product Datasheet

Human Sarcoplasmic Reticulum Histidine-rich Calcium-binding Protein ELISA Kit (orb1210149)

Description	Human Sarcoplasmic Reticulum Histidine-rich Calcium-binding Protein, HRC ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	HRC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	80,244
Uniprot ID	P23327
Solubility (25°C)	16ng/ml
Sensitivity	0.019ng/ml
Entrez	3270
Expiration Date	Please enquire.