

Product Datasheet

Human Guanylin ELISA Kit (orb1210152)

Description	Human Guanylin, GUCA2A ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	GUCA2A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	12,388
Uniprot ID	Q02747
Solubility (25°C)	12.8ng/ml
Sensitivity	0.017ng/ml
Entrez	2980
Expiration Date	Please enquire.