

Product Datasheet

Human Gastrokine 2 ELISA Kit (orb1210162)

Description	Human Gastrokine 2, GKN2 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	GKN2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,487
Uniprot ID	Q86XP6
Solubility (25°C)	4800ng/L
Sensitivity	10.27ng/L
Entrez	200504
Expiration Date	Please enquire.