

Product Datasheet

Human Exosome Complex Component Rrp46 ELISA Kit (orb1210168)

Description	Human Exosome Complex Component Rrp46, EXOSC5 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	EXOSC5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,249
Uniprot ID	Q9NQT4
Solubility (25°C)	2400ng/L
Sensitivity	5.87ng/L
Entrez	56915
Expiration Date	Please enquire.