

Product Datasheet

Human Epithelial Membrane Protein 3 ELISA Kit (orb1210177)

Description	Human Epithelial Membrane Protein 3, EMP3 ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	EMP3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,429
Uniprot ID	P54852
Solubility (25°C)	16ng/ml
Sensitivity	0.015ng/ml
Entrez	2014
Expiration Date	Please enquire.