

Product Datasheet

Human Homeobox Protein Engrailed-2 ELISA Kit (orb1210180)

Description	Human Homeobox Protein Engrailed-2, EN2 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	EN2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	34,211
Uniprot ID	P19622
Solubility (25°C)	48ng/ml
Sensitivity	0.048ng/ml
Entrez	2020
Expiration Date	Please enquire.