

Product Datasheet

Human Dynein Heavy Chain 11 ELISA Kit (orb1210183)

Description	Human Dynein Heavy Chain 11, Axonemal, DNAH11 ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	DNAH11
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	520,369
Uniprot ID	Q96DT5
Solubility (25°C)	32ng/ml
Sensitivity	0.023ng/ml
Entrez	8701
Expiration Date	Please enquire.